


## Product Datasheet

### CALM1 antibody (orb542937)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CALM1                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                           |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Tested Applications</b> | ICC, IF, IHC                                                |
| <b>Immunogen</b>           | A synthetic peptide of human Calmodulin 1/2/3(NP_008819.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.            |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                 |
| <b>Note</b>                | For research use only                                       |
| <b>Isotype</b>             | IgG                                                         |
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>MW</b>                  | 16kDa                                                       |
| <b>Uniprot ID</b>          | <b>P0DP24, P0DP25, P0DP23</b>                               |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:200, IF: 1:50-1:200                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                             |



Immunofluorescence analysis of HeLa cell...