

Product Datasheet

SNRNP70 antibody (orb542950)

Description	Rabbit polyclonal antibody to SNRNP70
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SNRNP70(NP_003080.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19kDa/39kDa/50kDa/51kDa
Uniprot ID	P08621
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.