

Product Datasheet

ZNF273 antibody (orb542962)

Description	Rabbit polyclonal antibody to ZNF273
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human ZNF273(NP_066971.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58kDa/64kDa
Uniprot ID	Q14593
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.