

Product Datasheet

ANKHD1 antibody (orb542970)

Description	Rabbit polyclonal antibody to ANKHD1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ANKHD1(NP_060217.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59kDa/63kDa/64kDa/269kDa/277kDa
Uniprot ID	Q8IWZ3
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.