

Product Datasheet

DDX27 antibody (orb542971)

Description	Rabbit polyclonal antibody to DDX27
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human DDX27(NP_060365.7).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89kDa
Uniprot ID	Q96GQ7
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.