

Product Datasheet

AMICA1 antibody (orb542982)

Description	Rabbit polyclonal antibody to AMICA1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human AMICA1(NP_694938.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa/40kDa/43kDa/44kDa
Uniprot ID	Q86YT9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.