

Product Datasheet

ADK antibody (orb542991)

Description	Rabbit polyclonal antibody to ADK
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ADK(NP_001114.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34kDa/36kDa/38kDa/40kDa
Uniprot ID	P55263
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.