

Product Datasheet

EGR2 antibody (orb542993)

Description	Rabbit polyclonal antibody to EGR2
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human EGR2(NP_000390.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa/50kDa
Uniprot ID	P11161
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.