

Product Datasheet

ADAMTS4 antibody (orb542999)

Description	Rabbit polyclonal antibody to ADAMTS4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein of human ADAMTS4 (NP_005090.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	075173
Entrez	9507
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.