

Product Datasheet

RABGAP1L antibody (orb543000)

Description	Rabbit polyclonal antibody to RABGAP1L
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RABGAP1L(NP_055672.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8kDa, 29kDa, 41-43kDa, 68kDa, 90-94kDa
Uniprot ID	B7ZAP0, Q5R372
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.