

Product Datasheet

NAP1L4 antibody (orb543017)

Description	Rabbit polyclonal antibody to NAP1L4
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NAP1L4 (NP_005960.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q99733
Entrez	4676
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.