



Product Datasheet

PFDN4 antibody (orb543021)


Description	Rabbit polyclonal antibody to PFDN4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PFDN4(NP_002614.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa
Uniprot ID	Q9NQP4
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human CO...



Immunohistochemical staining of rat live...