

Product Datasheet

SERPINB4 antibody (orb543022)

Description	Rabbit polyclonal antibody to SERPINB4
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human SERPINB4(NP_002965.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45KDa
Uniprot ID	P48594
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.