

Product Datasheet

SEC23IP antibody (orb543033)

Description	Rabbit polyclonal antibody to SEC23IP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SEC23IP(NP_009121.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102kDa/111kDa
Uniprot ID	Q9Y6Y8
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.