

Product Datasheet

FNDC3A antibody (orb543035)

Description	Rabbit polyclonal antibody to FNDC3A
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FNDC3A(NP_001265367.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125kDa/131kDa
Uniprot ID	Q9Y2H6
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.