

Product Datasheet

MPRIP antibody (orb543037)

Description	Rabbit polyclonal antibody to MPRIP
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MPRIP(NP_055949.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	114kDa/116kDa/118kDa
Uniprot ID	Q6WCQ1
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.