

Product Datasheet

TAF7L antibody (orb543046)

Description	Rabbit polyclonal antibody to TAF7L
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TAF7L (NP_079161.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q5H9L4
Entrez	54457
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.