

Product Datasheet

RC3H1 antibody (orb543064)

Description	Rabbit polyclonal antibody to RC3H1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RC3H1(NP_742068.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	126KDa
Uniprot ID	Q5TC82
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.