

Product Datasheet

ADGRA3 antibody (orb543066)

Description	Rabbit polyclonal antibody to ADGRA3
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ADGRA3(NP_660333.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	146kDa
Uniprot ID	Q8IWK6
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.