

Product Datasheet

CD163L1 antibody (orb543070)

Description	Rabbit polyclonal antibody to CD163L1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CD163L1(NP_777601.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79kDa/155kDa/159kDa/160kDa
Uniprot ID	Q9NR16
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.