

Product Datasheet

MIA3 antibody (orb543071)

Description	Rabbit polyclonal antibody to MIA3
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MIA3(NP_940953.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56kDa/87kDa/206kDa/213kDa
Uniprot ID	Q5JRA6
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.