

Product Datasheet

BMP2K antibody (orb543096)

Description	Rabbit polyclonal antibody to BMP2K
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BMP2K(NP_060063.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73kDa/75kDa/129kDa
Uniprot ID	Q9NSY1
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.