

Product Datasheet

NMRAL1 antibody (orb543097)

Description	Rabbit polyclonal antibody to NMRAL1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NMRAL1 (NP_065728.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9HBL8
Entrez	57407
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.