

Product Datasheet

NOD2 antibody (orb543102)

Description	Rabbit polyclonal antibody to NOD2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human NOD2(NP_071445.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100KDa
Uniprot ID	Q9HC29
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.