



Product Datasheet

IGHMBP2 antibody (orb543116)


Description	Rabbit polyclonal antibody to IGHMBP2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human IGHMBP2(NP_002171.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	109kDa
Uniprot ID	P38935
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human ma...



Immunohistochemical staining of human lu...