

Product Datasheet

NEK1 antibody (orb543117)

Description	Rabbit polyclonal antibody to NEK1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NEK1(NP_036356.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa/135kDa/138kDa/141kDa/142kDa/145kDa
Uniprot ID	Q96PY6
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.