

Product Datasheet

Human STK32A protein (orb54312)

Description	Recombinant human STK32A protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	N-terminal GST-tagged: N-terminal GST-tagged1-229AA: 1-166AAFull Length of BC016640 : Full Length of Isoform 2
Protein Sequence	MGANTSRKPPVFDENEDVNFDFEILRAIGKGSFGKVCIVQKNDTKKMYAMKYM NKQKCVERNEVRNVFKELQIMQGLEHPFLVNLWYSFQDEEDMFMVVDLLLGGD LRYHLQQNVHFKEETVKLFICELVMALDYLQNQRRIHRDMKPDNILLDEHDTWLS YKSH
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	46.8 kDa
Uniprot ID	Q8WU08
Expiration Date	6 months from date of receipt.