

Product Datasheet

PDE11A antibody (orb543128)

Description	Rabbit polyclonal antibody to PDE11A
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human PDE11A (NP_058649.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	105kDa
Uniprot ID	Q9HCR9
Entrez	50940
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.