

Product Datasheet

PDE11A antibody (orb543129)

Description	Rabbit polyclonal antibody to PDE11A
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PDE11A(NP_058649.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa/65kDa/78kDa/104kDa
Uniprot ID	Q9HCR9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.