

Product Datasheet

PADI4 antibody (orb543140)

Description	Rabbit polyclonal antibody to PADI4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PADI4(NP_036519.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	74kDa
Uniprot ID	Q9UM07
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.