

Product Datasheet

SV2B antibody (orb543166)

Description	Rabbit polyclonal antibody to SV2B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SV2B(NP_055663.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa/77kDa
Uniprot ID	Q7L1I2
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.