

Product Datasheet

Collagen IV antibody (orb543187)

Description	Rabbit Polyclonal antibody to Collagen IV
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human COL4A5(NP_203699.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	161kDa
Uniprot ID	P29400
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.