

Product Datasheet

GRINA antibody (orb543196)

Description	Rabbit polyclonal antibody to GRINA
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human GRINA(NP_000828.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53kDa
Uniprot ID	Q7Z429
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.