

Product Datasheet

PON3 antibody (orb543212)

Description	Rabbit polyclonal antibody to PON3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human PON3 (NP_000931.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	39kDa
Uniprot ID	Q15166
Entrez	5446
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.