

Product Datasheet

TAF6 antibody (orb543216)

Description	Rabbit polyclonal antibody to TAF6
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human TAF6(XP_006716164.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	71kDa/72kDa/76kDa
Uniprot ID	P49848
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.