

Product Datasheet

NR2F1 antibody (orb543217)

Description	Rabbit polyclonal antibody to NR2F1
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human NR2F1 (NP_005645.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P10589
Entrez	7025
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.