

## Product Datasheet

### PARP2 antibody (orb543230)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PARP2                 |
| <b>Reactivity</b>          | Mouse, Rat                                          |
| <b>Conjugation</b>         | Unconjugated                                        |
| <b>Tested Applications</b> | WB                                                  |
| <b>Immunogen</b>           | A synthetic peptide of human PARP2(NP_001036083.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.         |
| <b>Note</b>                | For research use only                               |
| <b>Isotype</b>             | IgG                                                 |
| <b>Clonality</b>           | Polyclonal                                          |
| <b>MW</b>                  | 64kDa/66kDa                                         |
| <b>Uniprot ID</b>          | <b>Q9UGN5</b>                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |