

Product Datasheet

SPINK5 antibody (orb543231)

Description	Rabbit polyclonal antibody to SPINK5
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SPINK5(NP_001121170.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	104kDa/120kDa/124kDa
Uniprot ID	Q9NQ38
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.