

Product Datasheet

KATNA1 antibody (orb543232)

Description	Rabbit polyclonal antibody to KATNA1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human KATNA1(NP_001191005.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35kDa/55kDa
Uniprot ID	075449
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.