

Product Datasheet

SOX21 antibody (orb543233)

Description	Rabbit polyclonal antibody to SOX21
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SOX21(NP_009015.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	Q9Y651
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.