

## Product Datasheet

### HDAC7 antibody (orb543248)

|                            |                                                  |
|----------------------------|--------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to HDAC7              |
| <b>Reactivity</b>          | Human, Rat                                       |
| <b>Conjugation</b>         | Unconjugated                                     |
| <b>Tested Applications</b> | ICC, IF, IHC                                     |
| <b>Immunogen</b>           | A synthetic peptide of human HDAC7(NP_056216.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.      |
| <b>Note</b>                | For research use only                            |
| <b>Isotype</b>             | IgG                                              |
| <b>Clonality</b>           | Polyclonal                                       |
| <b>MW</b>                  | 46-66kDa, 99-108kDa                              |
| <b>Uniprot ID</b>          | <b>Q8WUI4</b>                                    |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:200, IF: 1:50-1:200                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                  |