

Product Datasheet

SOX18 antibody (orb543249)

Description	Rabbit polyclonal antibody to SOX18
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human SOX18(NP_060889.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	P35713
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.