

Product Datasheet

ZCCHC10 antibody (orb543251)

Description	Rabbit polyclonal antibody to ZCCHC10
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZCCHC10(NP_001287745.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	Q8TBK6
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.