

Product Datasheet

YY2 antibody (orb543273)

Description	Rabbit polyclonal antibody to YY2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human YY2(NP_996806.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	O15391
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.