

Product Datasheet

TERT antibody (orb543274)

Description	Rabbit polyclonal antibody to TERT
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human TERT(NP_937983.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	88kDa/90kDa/120kDa/126kDa
Uniprot ID	O14746
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.