

Product Datasheet

TTN antibody (orb543296)

Description	Rabbit polyclonal antibody to TTN
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human TTN(NP_596870.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	3816KDa
Uniprot ID	Q8WZ42
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.