

Product Datasheet

SNIP1 antibody (orb543305)

Description	Rabbit polyclonal antibody to SNIP1
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human SNIP1(NP_078976.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46kDa
Uniprot ID	Q8TAD8
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.