



Product Datasheet

CHRM1 antibody (orb543358)


Description	Rabbit polyclonal antibody to CHRM1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CHRM1(NP_000729.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51KDa
Uniprot ID	P11229
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of OVCAR-3
C...



Immunofluorescence analysis of NIH-3T3
C...



Immunofluorescence analysis of C6 cells
...