

Product Datasheet

NR4A2 antibody (orb543669)

Description	Rabbit polyclonal antibody to NR4A2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic Peptide of human NR4A2.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	67KDa
Uniprot ID	P43354
Entrez	4929
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.