

Product Datasheet

MAP1S antibody (orb543715)

Description	Rabbit polyclonal antibody to MAP1S
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MAP1S(NP_060644.4).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	109kDa/112kDa
Uniprot ID	Q66K74
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.